



182x208x46mm



- Main device with replaceable earbuds (M size)
- 3 x Additional pair of L, S & Mini size earbuds
- Cleaning Brush
- USB Cable

Item No.	GoodHearing A1
Max. Output @ SPL 90 dB	135 ± 5 dB
Max Gain	60 ± 5 dB
Continuous Use	Up to 10 hours per charge
Listening Mode	3 Modes
Suitable For	People with moderate to severe hearing loss
Storage Temperature	14°F ~ 140°F
Storage Humidity	30% ~ 85%
Dimension	160 x 190 x 11 mm
Weight	35g
Place of Origin	Taiwan



Copyright © Digisine 2020.V.1.0



*Nguyễn Hoàng Phương*

Mimitakara®

Rechargeable Hearing Aid

Digital Active Type

Mimitakara®

Digital Active Type

## Rechargeable Hearing Aid

Dual microphone  
Patented audio technology  
Smart control with Bluetooth & App



Digital 8 Channel  
Suitable for both ear

OPEN



ISO 13485 CE 2460

FDA Registration# 3005053820  
Owner/Operator# 9071124

Digital Active Type

Rechargeable Hearing Aid

Mimitakara®