

135x50x180mm



Mimitakara@myHearing I1 (ITC Type) Technical Data

Type	ITC (I)	
	2 ccm coupler	Ear simulator
Output sound pressure level at 1.6 kHz	-	121 dB SPL
Peak	116 dB SPL	127 dB SPL
HFA-OSPL 90	102 dB SPL	-
Gain		
Full on gain (FOG) at 1.6 kHz	-	55 dB
Full on gain (Peak)	55 dB	65 dB
HFA-FOG	49 dB	-
Reference test gain	35 dB	45 dB
Frequency, noise and directivity		
Frequency range	100-8100 Hz	120-8100 Hz
Equivalent input noise	21 dB SPL	22 dB SPL
Total harmonic distortion at 500 / 800 / 1600 Hz	2 / 2 / 1 %	3 / 4 / 5 %
Battery		
Battery voltage	1.3 V	
Battery current drain	1.0 mA	
Battery life (cell zinc air)	120h	

I1 : Right ear



4 719878 170208

I1 : Left ear



4 719878 170215

40X30 白底黑字貼紙



Nguyễn Hoàng Phương



局部凸

