

Label Set 08717010190 ACTB12 - Elecsys Active B12 Version 2
 eCDS Label Set Item No. 1485 V 1
 Release Date 12/1/2021 12:34:25 PM
 Special instructions:

Elecsys Active B12

cobas e analyzers

REF 08717010190

COBAS, COBAS E and ELECSYS are trademarks of Roche.

Roche Diagnostics GmbH
 Sandhofer Str. 116
 D-68305 Mannheim

07020299001(1) Distribution in USA by:
 Roche Diagnostics,
 Indianapolis, IN
 Made in Germany

Roche

H317
 P261,P272,P280
 P333+P313,P362+P364
 P501

CONTENT 100

M	6.5 mL
R1	9 mL
R2	9 mL

IVD 0123 2-8 °C

dialog.roche.com

01 **cobas**[®]

Pos.1 - 09221026700 - V1

ACTB12

REF (240)08717010190

GTIN
 (01)07613336163429

UDI

LOT (10)00000000

2000-01-31

2000-01-01

01 **Roche**

Pos.2 - 09221034700 - V1


Label Set 08717010190 ACTB12 - Elecsys Active B12 Version 2
eCDS Label Set Item No. 1485 V 1
Release Date 12/1/2021 12:34:25 PM
Special instructions:


Elecsys Active B12

REF 08717010190


For **USA**: **CONTENT**

- M** Streptavidin-coated microparticles 0.72 mg/mL
- R1** Biotinylated monoclonal anti-holotranscobalamin antibody (mouse) 2 mg/L
- R2** Monoclonal anti-transcobalamin antibody (mouse) labeled with ruthenium complex 0.5 mg/L



Rx only 01 

Pos.3 - 09221085700 - V1



ACTB12


Elecsys **1550**


08717010


- M** 6.5 mL
- R1** 9 mL
- R2** 9 mL


1

LOT 00000000

 2000-01-31

 2-8 °C

IVD 

0101 

000004294667296

Pos.4 - 09221042700 - V1

Label Set 08717010190 ACTB12 - Elecsys Active B12 Version 2

eCDS Label Set Item No. 1485 V 1

Release Date 12/1/2021 12:34:25 PM

Special instructions:

ACTB12

167

1

LOT 00000000

2000-01-31

H317

000004294967296

0101

Roche

07021127(1)

The label features a large vertical barcode on the left side. To the right of the barcode, the text 'ACTB12' is prominently displayed at the top. Below it, the number '167' is enclosed in a small box. Further down, the number '1' is printed. The 'LOT' number '00000000' and the expiration date '2000-01-31' are shown with their respective icons. A red diamond warning symbol with an exclamation mark is present, with the code 'H317' to its left and '07021127(1)' to its right. A vertical alphanumeric string '000004294967296' is located to the left of the warning symbol. At the bottom left, the number '0101' is printed, and the Roche logo is at the bottom right. A small square QR code is positioned above the Roche logo.

Pos.5 - 09221069700 - V1

Nhãn sản phẩm
Elecsys Active B12 (08717010190)

Thuốc thử xét nghiệm **định lượng** vitamin B12
hoạt tính (holotranscobalamin)

Elecsys Active B12

Mã sản phẩm: 08717010190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:

Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , Số lô: xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991