



■ 燙銀

96x140x32mm



Product Illustration:



- ① Hearing Device x 1
  - ② Portable Case x 1
  - ③ Cleaning Brush x 1
  - ④ Ear Cap (S,M,L x 1)
  - ⑤ Trial Battery x 2
- ※ Included Zn Air Battery 312(PRA1) is for first time testing use.

How to Use the Hearing Device?

1. Open the portable case and take out the hearing device. Insert the battery into the battery compartment. Press down on the battery cap.
2. The hearing device comes with 3 different ear cap sizes (large/medium /small). Please select the size that best fits your ear.
3. Use your finger to grasp the microphone and block the hole. Place the hearing device into your ear.
4. Adjust the volume to a suitable level.  
If you hear a harsh sound, it means the volume setting is too high or the hearing device is not placed in your ear properly. Adjust the position of the hearing device or reduce the volume setting slightly to resolve this problem.
5. Take the hearing device out of your ear when you finish using it and place it in the portable case, close the case will automatically keep the battery isolated.



Made in Taiwan



4 711914 557612

Mimitakara® JAPAN  
Digital Hearing Device  
ITC Type Battery Operated

Mimitakara® JAPAN  
Digital Hearing Device  
ITC Type Battery Operated

ISO 13485 CE 2460  
FDA Registration# 3005053820  
Owner/Operator# 9071124



Mimitakara® JAPAN  
Digital Hearing Device  
ITC Type Battery Operated



*Nguyễn Hoàng Phương*



# 96x140x16.5mm

灰色代表燙銀

## Mimitakara

Mimitakara most respected hearing aid and hearing amplifier brand in Taiwan, Japan, Europe & USA. Mimitakara believes in QUALITY and SERVICE as the company's core foundation. We are the World's First Cloud Hearing Service and are a pioneer in the segment. We design, develop and manufacture all our products ourselves alongside our trusted partners in Taiwan, to ensure that our clients and their family to have the best hearing experience. Mimitakara believes EVERYONE DESERVES TO HEAR and brings high-tech hearing aids and amplifiers to you at factory cost by eliminating intermediary markups with international standards and certification.

FDA Registration 3005053820

JFDA BG10600253

ISO 13485

CE 2460



Copyright © Digisine 2022.V.1.6



4 719878 171403

Mimitakara® ITC Type Hearing Device

## Mimitakara®

### ITC Type

## Hearing Device

ISO 13485 CE 2460

FDA Registration# 3005053820  
Owner/Operator# 9071124



Mimitakara® ITC Type Hearing Device

