



## Document Information:

Document Author: Limburger, Marion {DQXP~Penzberg}  
Business Area / Unit: Roche Professional Diagnostics  
Confidentiality: Confidential  
Document Class: Device Master Record  
Document Type: Device Master Record  
Document Creator: Limburger, Marion {DQXP~Penzberg}  
Document Lifecycle Status: Signed  
Valid From: 22-Jan-2021 10:24:49 (UTC)  
Valid To:  
Document Title: C\_08\_RSP\_07026684190\_EN  
Document Number: 0000000000001200000173293  
Document Version: 03  
Template: No

Global Group: Regulatory Affairs  
Global SubGroup: Specification Prod. (Registr.)  
Local Group: Reg. Specification  
Language: English  
Site: RDG Germany  
Department: C\_DOMS  
Document Applies To: C\_Standardization/ QC Penzberg  
Document Description: Elecsys ACTH - "Only for use in Regulatory Affairs"

## Electronic Signatures:

Signed By: limburgm (Marion Limburger {DQXPFB})  
Role: Author  
Signature Differentiation: Quality Assurance  
Signed Date: 19-Jan-2021 06:53:30 (UTC)

Signed By: wickerts (Sabine Wickert {DOMSW})  
Role: Reviewer  
Signed Date: 21-Jan-2021 08:46:25 (UTC)

Signed By: tassionir (Richard Tassoni {DOMS})  
Role: Approver  
Signed Date: 22-Jan-2021 10:14:37 (UTC)

# SPECIFICATION

for

**Elecsys ACTH**

---

**0702 6684 190**

**0894 6728 190**

**0325 5751 190**

**0894 6710 190**

## Part 1 (M)

Appearance	sediment with clear supernatant
Streptavidin-coated microparticles	648 - 792 µg/mL

## Part 2 (R1)

Appearance	clear solution
Biotinylated monoclonal anti-ACTH antibody	proved by performance test
pH value	6.0 - 6.4

## Part 3 (R2)

Appearance	clear solution
Monoclonal anti-ACTH antibody labeled with ruthenium complex	proved by performance test
pH value	6.0 - 6.4

## Performance test

### Recovery of ACTH in

(based on assigned value)

PreciControl Multimarker Level 1	80 - 120 %
PreciControl Multimarker Level 2	80 - 120 %