



Document Information:

Document Author:	Schwaiger, Nicole {DOMS~Penzberg}
Business Area / Unit:	Roche Professional Diagnostics
Confidentiality:	Confidential
Document Class:	Device Master Record
Document Type:	Device Master Record
Document Creator:	Schwaiger, Nicole {DOMS~Penzberg}
Document Lifecycle Status:	Signed
Valid From:	09-May-2022 06:07:17 (UTC)
Valid To:	
Document Title:	C_08_RSP_08846715190_EN
Document Number:	0000000000001200000452860
Document Version:	02
Template:	No
Global Group:	Regulatory Affairs
Global SubGroup:	Specification Prod. (Registr.)
Local Group:	Reg. Specification
Language:	English
Site:	RDG Germany
Department:	C_DOMS
Document Applies To:	C_Standardization/ QC Penzberg
Document Description:	Elecsys Phospho-Tau (181P) CSF - "Only for use in Regulatory Affairs"

Electronic Signatures:

Signed By:	gehrker (Ralph Gehrke {DOMS})
Role:	Approver
Signed Date:	09-May-2022 06:01:40 (UTC)

SPECIFICATION

for

Elecsys Phospho-Tau (181P) CSF

0884 6715 190

0884 6693 190

Part 1 (M)

Appearance	sediment with clear supernatant
Streptavidin coated micro particles	648 - 792 µg/ml

Part 2 (R1)

Appearance:	clear solution
Biotinylated monoclonal anti-pTau antibody	proved by performance test
pH value	6.9 – 7.5

Part 3 (R2)

Appearance:	clear solution
Monoclonal anti-Tau antibody, labeled with ruthenium complex	proved by performance test
pH value	6.9 – 7.5

Performance test

Recovery of pTau (based on assigned value)

PreciControl 1 pTau	86 – 114 %
PreciControl 2 pTau	86 – 114 %