

Elecsys Phospho-Tau (181P) CSF

cobas e analyzers

H317
P261,P272,P280
P333+P313
P362+P364,P501



REF 08846693190

CONTENT Σ 60

COBAS, COBAS E
and ELECSYS are
trademarks of Roche

M 6.5 mL
R1 6.5 mL
R2 6.5 mL



 Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

070202959001 (1)
Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN
Made in Germany

IVD **CE** 0123  2-8 °C

01  dialog.roche.com



cobas[®]

Pos.1 - 09226451700 - V1

pTau

REF (240)08846693190


GTIN

(01)07613336167656

UDI



LOT (10)0000000

 2000-01-31

 2000-01-01

02



Pos.2 - 09226460700 - V2

Elecsys Phospho-Tau (181P) CSF

REF 08846693190

For USA: CONTENT

- M Streptavidin-coated microparticles 0.72 mg/mL
- R1 Biotinylated monoclonal anti-pTau antibody (rabbit/mouse) 2.5 mg/L
- R2 Monoclonal anti-Tau antibody (mouse) labeled with ruthenium complex 2.0 mg/L



Rx only

02



Pos.3 - 09226478700 - V2

pTau
Elecsys 1670
08846693

M	6.5 mL
R1	6.5 mL
R2	6.5 mL

1
LOT 00000000
2000-01-31
2-8 °C

IVD

0101


000004294967296

Pos.4 - 09226494700 - V1


pTau 179

1


LOT 00000000


 2000-01-31

H317

 07021127(1)

00000429467296



0101 

Pos.5 - 09226486700 - V1

Nhãn sản phẩm
Elecsys Phospho-Tau (181P) CSF (08846693190)

Thuốc thử xét nghiệm định lượng protein Tau
phosphoryl hoá
Elecsys Phospho-Tau (181P) CSF

Mã sản phẩm: 08846693190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:

Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem , Số lô: xem **LOT**

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991