

TEM NHÃN HIỆU HÀNG HÓA

NHÃN HIỆU: HUANG YU

NHÃN HIỆU: ZMC

ỐNG LẤY MÁU XÉT NGHIỆM EDTA K3

 HUANG YU	VACUUM TUBE EDTA K3 TUBE 2ml
Color Code: Purple	Mfg Date: 20170402
Size: 13x75mm	LOT: 20170402
QTY: 100pcs	⏰: 20190401
 Made In China	 6 901046 601456

730

size: 110x60mm

ỐNG LẤY MÁU XÉT NGHIỆM EDTA K3



 ZMC	VACUUM TUBE EDTA K3 TUBE 2ml
Color Code: Purple	Mfg Date: 20170511
Size: 13x75mm	LOT: 20170511
QTY: 100pcs	⏰: 20190510
 Made In China	 6 901046 601456

**ỐNG LẤY MÁU XÉT NGHIỆM
SODIUM CITRAT 3,8 % (1:9)**

	VACUUM TUBE 3.8% Sodium Citrate 2ml
Color Code: Green	Mfg Date: 20170402
Size: 13x75mm	LOT: 20170402
QTY: 100pcs	⌚ : 20190401
 Made In China	

730


**ỐNG LẤY MÁU XÉT NGHIỆM
SODIUM CITRAT 3,8 % (1:9)**

	VACUUM TUBE PT Tube 3.8% Sodium Citrate 1:9 2ml
Color Code: Green	Mfg Date: 20170126
Size: 13x75mm	LOT: 20170126
QTY: 100pcs	⌚ : 20190125
 Made In China	

ỐNG LẤY MÁU XÉT NGHIỆM HEPARIN LITHIUM

	VACUUM TUBE Heparin Lithium 2ml
Color Code: Black	Mfg Date: 20170402
Size: 13x75mm	LOT: 20170402
QTY: 100pcs	⌚ : 20190401
 Made In China	

730




ỐNG LẤY MÁU XÉT NGHIỆM HEPARIN LITHIUM

	VACUUM TUBE Heparin Lithium 2ml
Color Code: Black	Mfg Date: 20170511
Size: 13x75mm	LOT: 20170511
QTY: 100pcs	⌚ : 20190510
 Made In China	

ĐẦU CÔN XANH




ĐẦU CÔN XANH



ĐẦU CÔN VÀNG



ĐẦU CÔN VÀNG



ỐNG ĐỰNG MẪU XÉT NGHIỆM 13X75



ỐNG ĐỰNG MẪU XÉT NGHIỆM 13X75



QUE TẮM BÔNG VÔ TRÙNG:



LAM KÍNH 7102:



LAM KÍNH 7105:



TÚI NƯỚC TIỂU :



DÂY OXY 2 NHÁNH:

▲
TEAR HERE

HƯƠNG VU
Nasal Oxygen Cannula

With Chamber and Tubing   0123

Directions For Use:

1. Attach the oxygen supply tubing to the gas source and set the oxygen to the prescribed flow.
2. Check for gas flow throughout the device.
3. Insert the nasal tips into the nostrils with the two plastic tubes over the ears and under the chin.
4. Gently adjust the plastic side until the cannula is secure.

NOTE: Nasal tips may be trimmed with scissors to fit smaller patients.

Caution:


- For single use. Discard after use.
- Do not use if package is open or damaged.
- Do not store at extreme temperature or humidity. Store in dry place and without sunlight.
- Ensure that all connections are secure and air circulates freely through the tubing.

SIZE: XS S M L

170329 2017-09 2022-08

IMPORTER: HOANG QUAN MEDICAL EQUIPMENT JOINT STOCK COMPANY
Address: 181 Kim Hoa st., Phung Liet w., Dong Da dist., Hanoi, Vietnam
Tel: 843 972 1049/9912 85 1186
MANUFACTURER: ZHEJIANG MEDICINES & HEALTH PRODUCTS YE CO., LTD

KIM CHÍCH MÁU:



KIM BƯỚM CHÂN KHÔNG:

↑ LOT NO./MFG.DATE

Disposable Butterfly Blood Collection Needle

HUANG YU  **SIZE: 21Gx1"**
Sterile, Non-toxic, Non-pyrogenic
Do not use if the package is open or damaged
For single use only, discard after use
Valid for 2 years

STERILE EO 


DIRECTION FOR USE:

1. Open package and remove needle set, then remove adapter needle protector
2. Prep patient
3. Remove needle protector and perform venipuncture
4. Attach vacuum tube
5. Collect specimen

↓ OPEN HERE 

↑ LOT NO./MFG.DATE

Disposable Butterfly Blood Collection Needle

HUANG YU  **SIZE: 22Gx1"**
Sterile, Non-toxic, Non-pyrogenic
Do not use if the package is open or damaged
For single use only, discard after use
Valid for 2 years

STERILE EO 

DIRECTION FOR USE:

1. Open package and remove needle set, then remove adapter needle protector
2. Prep patient
3. Remove needle protector and perform venipuncture
4. Attach vacuum tube
5. Collect specimen

↓ OPEN HERE 

Disposable Butterfly Blood Collection Needle

SIZE: 23G

100Pcs/Bag

STERILE EO

CE

