

## INSTRUCTIONS FOR USE

Equipment Name : Underarm crutch



DY05925L (L); DY05925L (M); DY05925L (S); DY05925L (XS); DY05925L (2) (M);DY05925L (3) (S); DY05925L (3) (M); DY 05925L (3) (L); DY05925L (4) (L);DY05925L (4) (M); DY05925L (4) (S); DY05925L (5) (S); DY05925L (5) (M); DY05925L (5) (L); DY05925L (6) (S); DY 05925L (6) (M); DY05925L (6) (L), DY05935 (L); DY05935 (M; DY05935 (S); DY05935 (2) (M); DY05943L

Equipment Name : Forearm crutch



DY 05933L; DY05933L (S); DY05933L (2); DY 05933L (3); DY05933L (4); DY059331L;  
DY059331L (S); DY059331L (2); DY059331LB; DY059331L (3); DY05937L;  
DY05937LB; DY05937LN; DY05937L (3); DY05923L; DY05923L (2)

Used in clinics, medical examination and treatment facilities, in families of disabled people, elderly people with functional injuries in their legs, people with reduced mobility.

B1: Adjust the crutch height to suit the patient's height and comfortable grip

B2: Check that the pins are firmly locked to keep the length of the crutches when in use  
The circulation registration establishment commits that the above-declared contents are true and takes responsibility before law for the above-declared information.