

Label Set 07026773190 AHAVIGM - Elecsys Anti-HAV IgM Version 5
 eCDS Label Set Item No. 1109 V 1
 Release Date 9/27/2021 5:50:22 AM
 Special instructions:

Elecsys Anti-HAV IgM

cobas e analyzers

REF 07026773190

GTIN 04015630939855

COBAS, COBAS E and
ELECSYS are trademarks
of Roche.

Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

07026773190(1)
Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN
Made in Germany

CONTENT 300

R1	21.0 mL	H317 P261 P272 P280 P333+P313 P362+P364 P501
R2	21.0 mL	
M	12.4 mL	
Cal1	0.67 mL	
Cal2	0.67 mL	

LOT 00000000

2000-01-31

IVD 0123 2-8 °C

03 dialog.roche.com

cobas[®]

Pos.1 - 07617763700 - V3

AHAVIGM

REF (240)07026773190

GTIN
(01)04015630939855

UDI

LOT (10)00000000

2000-01-31

2000-01-01

04

Pos.2 - 07617771700 - V4

Label Set 07026773190 AHAVIGM - Elecsys Anti-HAV IgM Version 5
eCDS Label Set Item No. 1109 V 1
Release Date 9/27/2021 5:50:22 AM
Special instructions:

Elecsys Anti-HAV IgM

REF 07026773190

For **USA**: **CONTENT**


R1 Monoclonal anti-HAV antibody (mouse)
labeled with ruthenium complex 0.15 µg/mL;
anti-human-Fdy antibody (sheep) 0.04 mg/mL


R2 Biotinylated monoclonal anti-h-IgM antibody
(mouse) 0.4 µg/mL; HAV antigen (cell culture)
25 U/mL (Roche units)

M Streptavidin-coated microparticles 0.72 mg/mL

Cal1 Human serum, negative for anti-HAV IgM

Cal2 Anti-HAV IgM (human) approximately 5 U/mL
(Roche units) in human serum

Rx only 001 



Pos.3 - 07617780700 - V1

AHAVIGM
07026773190

R1
R2
M




IVD
LOT
00000000


2000-01-31


01

Pos.4 - 07617879700 - V1


Label Set 07026773190 AHAVIGM - Elecsys Anti-HAV IgM Version 5
eCDS Label Set Item No. 1109 V 1
Release Date 9/27/2021 5:50:22 AM
Special instructions:




Pos.5 - 07617887700 - V3



Pos.6 - 07617798700 - V1




Pos.7 - 07617836700 - V1



Pos.8 - 07617810700 - V1

Label Set 07026773190 AHAVIGM - Elecsys Anti-HAV IgM Version 5
eCDS Label Set Item No. 1109 V 1
Release Date 9/27/2021 5:50:22 AM
Special instructions:



Pos.9 - 07617844700 - V1

Nhãn sản phẩm
Elecsys Anti-HAV IgM (07026773190)

Thuốc thử xét nghiệm định tính kháng thể IgM
kháng HAV
Elecsys Anti-HAV IgM

Mã sản phẩm: 07026773190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:

Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem , Số lô: xem **LOT**

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991