



VENTANA®

*a member of the Roche group*

## QUALITY CONTROL SPECIFICATION

---

**Product name: Diastase Kit**

**Ventana REF : 860-004**

**Roche Number: 05279208001**

Functional Staining on  
liver tissue

Pas-Hematoxylin stains glycogen bright magenta

Diastase slides should have a trace to no magenta  
present

Negative control – Cells where the cell nuclei must  
stain purple

 Vincent A. Faddis 22-NOV-2011  
**Quality Department Manager** **Date**