

UREAL

cobas c 111

REF 04657616190

Urea/BUN

COBAS and COBAS C are trademarks of Roche.

Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim


Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN


Made in Germany

02


CONTENT
R1 **REAGENT** 4 x 7.0 mL

Σ 4 x 100

IVD **CE** 0123  2-8 °C

 dialog.roche.com

Roche **cobas**





UREAL


REF (240)04657616190

GTIN (01)04015630921966

LOT (10)00000000

UDI 

 2000-01-31

 2000-01-01

04 **Roche**

For **USA**:


REF 04657616190

CONTENT

R1: TRIS buffer 220 mmol/L, pH 8.6; 2-Oxoglutarate 73 mmol/L; NADH 2.5 mmol/L; ADP 6.5 mmol/L; Urease (jack bean) ≥ 300 μkat/L; GLDH (bovine liver) ≥ 80 μkat/L

Rx only

02 **Roche**




1 **UREAL** **1**

04657616


7.0 mL


REAGENT

 2-8 °C


IVD

LOT
00000000

 2000-01-31



002



12311294567812120000001

**Nhãn sản phẩm
UREAL (04657616190)**

**Thuốc thử xét nghiệm định lượng urea/urea
nitrogen
UREAL**

Mã sản phẩm: 04657616190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:

Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , **Số lô:** xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991