



REF 09P5220

KXXXXV
GXXXX/RXX



Thuốc thử xét nghiệm định tính và/hoặc bán định lượng benzodiazepines

Alinity c Benzodiazepines Reagent Kit

Số lưu hành: XXXX-XXXX-XXXX

Chủ sở hữu số lưu hành: Văn phòng Đại diện Abbott Laboratories GmbH tại Hà Nội. Địa chỉ: tầng 7&8, tháp A, toà nhà Handi Resco, số 521 Kim Mã, Ba Đình, Hà Nội, Việt Nam.

Chủ sở hữu trang thiết bị: Microgenics Corporation, 46500, Kato Road Fremont, CA 94538, USA.

Xuất xứ: Hoa Kỳ.

Số lô và Ngày hết hạn: Xem "**LOT**" và "**Exp.**" trên nhãn gốc.


Nhiệt độ bảo quản: 2-8°C

Thông tin cảnh báo, Hướng dẫn sử dụng và các thông tin

khác: Xem Hướng dẫn sử dụng kèm theo hộp hoặc tài trên trang www.corelaboratory.abbott/IFU

Benz urine

Alinity c Benzodiazepines Reagent Kit



Benz urine

REF 09P5220

Σ 2 x 250



www.abbottdiagnostics.com/IFU




LOT

R1 2 x 32.2 mL

R2 2 x 31.8 mL

CONTAINS: AZIDE

G99829R01



Microgenics Corporation
4650 Kato Road
Fremont, CA 94538 USA

EC REP

B-R-A-H-M-S GmbH
Neuendorfstrasse 25
16761 Hennigsdorf
Germany

DISTRIBUTED BY

Abbott Laboratories

IVD

PRODUCT OF USA



INFORMATION FOR USA ONLY

Rx ONLY

The Alinity c Benzodiazepines assay is a homogeneous enzyme immunoassay intended for the qualitative and/or semiquantitative determination of the presence of benzodiazepines and their metabolites in human urine at a cutoff concentration of 200 ng/mL (0.700 μ mol/L) on the Alinity c analyzer.

- R1** Active ingredients: Anti-benzodiazepine polyclonal antibodies (sheep) (\leq 3 mg/mL), G6P (\leq 0.0062 g/mL), NAD (\leq 0.0076 g/mL). Inactive ingredients: TRIS buffer, bovine serum albumin (BSA). Preservative: sodium azide (\leq 0.09%).
- R2** Active ingredient: Benzodiazepine derivative labeled with G6PDH (\leq 0.8 mg/mL). Inactive ingredients: TRIS buffer, BSA. Preservative: sodium azide (\leq 0.09%).

Distributed in the USA by Abbott Laboratories, Abbott Park, IL 60064 USA

File: 09P52-TF-20_Eng_POTe.indd
Template: AK4498Rgt-CC-01_Eng_POTe.indt
Size: 181 mm x 58 mm
Stamp field: Refer to drawing
Colors: Text, Symbols: Black; Branding: Pantone 233c, Text reversed
Material: Refer to drawing

Drawing: AK4498