



PINMED

NINGBO PINMED INSTRUMENTS CO., LTD.

Use manual

Name: medical cotton wool ball

Specifications: 25g, 50g, 100g, 250g, 400g, 454g, 500g, 1000g

Caution: do not contact the wounded skin. Do not use the products when they are found to be yellow, mouldy, wet, and polluted.

Use directions :

The product is non-sterile. The users can dip in little amount of iodophor or alcohol before using the product when it's necessary.

Intended use: for medical dressing

Storage: store the product in cool and clean area, where is well-ventilated, without corrosive gas , with the relevant humidity of no more than 80%, temperature less than 50°C