



Document Information:

Document Author: Limburger, Marion {DQXP~Penzberg}
Business Area / Unit: Roche Professional Diagnostics
Confidentiality: Confidential
Document Class: Device Master Record
Document Type: Device Master Record
Document Creator: Limburger, Marion {DQXP~Penzberg}
Document Lifecycle Status: Effective
Valid From: 06-Dec-2019 11:44:09 (UTC)
Valid To:
Document Title: C_08_RSP_07027532190_EN
Document Number: 0000000000001200000195749
Document Version: 02
Template: No

Global Group: Regulatory Affairs
Global SubGroup: Specification Prod. (Registr.)
Local Group: Reg. Specification
Language: English
Site: RDG Germany
Department: C_DOMS
Document Applies To: C_Standardization/ QC Penzberg
Document Description: Elecsys IL-6 - "Only for use in Regulatory Affairs"

Electronic Signatures:

Signed By: limburgm (Marion Limburger {DQXPFB})
Role: Author
Signature Differentiation: Quality Assurance
Signed Date: 05-Dec-2019 08:13:21 (UTC)

Signed By: schrankr (Ralf Schrank {DOMSW})
Role: Reviewer
Signature Differentiation: Quality Control
Signed Date: 06-Dec-2019 11:40:31 (UTC)

Signed By: tassonir (Richard Tassoni {DOMS})
Role: Approver
Signed Date: 06-Dec-2019 11:41:30 (UTC)

SPECIFICATION

for

Elecsys IL-6

0702 7532 190

0901 5612 190

0510 9442 190

0510 9442 191

0901 5604 190

Part 1 (M)

Appearance

sediment with clear supernatant

Streptavidin-coated microparticles

648 - 792 µg/ml

Part 2 (R1)

Appearance

clear solution

Biotinylated monoclonal anti-IL-6 antibody

proved by performance test

pH-value

7.1 – 7.5

Part 3 (R2)

Appearance

clear solution

Monoclonal anti-IL-6 antibody ~ ruthenium complex

proved by performance test

pH-value

7.1 – 7.5

Performance test

Recovery of IL-6 in

(based on assigned value)

Control 1

PreciControl Multimarker Level 1

79 – 121 %

PreciControl Multimarker Level 2

79 – 121 %
