



Document Information:

Document Author: Limburger, Marion {DOMS~Penzberg}
Business Area / Unit: Roche Professional Diagnostics
Confidentiality: Confidential
Document Class: Device Master Record
Document Type: Device Master Record
Document Creator: Limburger, Marion {DOMS~Penzberg}
Document Lifecycle Status: Signed
Valid From: 13-Sep-2021 07:19:50 (UTC)
Valid To:
Document Title: C_08_RSP_05109469190_EN
Document Number: 0000000000001008018000860
Document Version: 03
Template: No

Global Group: Regulatory Affairs
Global SubGroup: Specification Prod. (Registr.)
Local Group: Reg. Specification
Language: English
Site: RDG Germany
Department: C_DOMS
Document Applies To: C_Standardization/ QC Penzberg
Document Description: IL-6 CalSet - "Only for use in Regulatory Affairs"

Electronic Signatures:

Signed By: limburgm (Marion Limburger {DOMSO})
Role: Author
Signature Differentiation: Quality Control
Signed Date: 06-Sep-2021 11:49:52 (UTC)

Signed By: wunderm6 (Michael Wunderlich {DOMSW})
Role: Reviewer
Signature Differentiation: Quality Control
Signed Date: 09-Sep-2021 16:19:28 (UTC)

Signed By: gehrker (Ralph Gehrke {DOMS})
Role: Approver
Signed Date: 13-Sep-2021 07:11:54 (UTC)

SPECIFICATION

for

IL-6 CalSet

0510 9469 190

Bottle 1 (IL-6 Cal1)

Appearance
IL-6 (human, recombinant)

lyophilized substance
approx. 18 pg/mL

Bottle 2 (IL-6 Cal 2)

Appearance
IL-6 (human, recombinant)

lyophilized substance
approx. 700 pg/mL

cobas e 801 analyzer: The exact lot-specific calibrator values are encoded in the electronic barcode and available via the cobas link.

All other analyzers: The exact lot-specific calibrator values are encoded in the barcode as well as printed on the enclosed (or electronically available) calibrator barcode sheet.