

IL-6 CalSet

Elecsys and
cobas e analyzers

REF 05109469190

CONTENT

Cal1 → 2 x 2.0 mL

Cal2 → 2 x 2.0 mL



COBAS, COBAS E
and ELECSYS are
trademarks of Roche.

Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

0644467001(1)
Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN
Made in Germany

IVD **CE** 0123 2-8 °C

dialog.roche.com

03



cobas[®]

Pos.1 - 07598874700 - V3

IL-6 CalSet

For USA: **CONTENT**
For Emergency Use Authorization only
Refer to Reagent Bulletin TP-00936

Cal1 IL-6 in equine serum matrix
Cal2 IL-6 in equine serum matrix

REF (240)05109469190

LOT (10)00000000

2000-01-31

2000-01-01

GTIN (01)04015630924233



UDI

Rx only

05



Pos.2 - 07598882700 - V5

IL6 Cal1
Elecsys → 2.0 mL
05109469

IVD 002 2-8 °C

LOT 00000000 2000-01-31

Pos.3 - 07598904700 - V2

IL6 Cal2
Elecsys → 2.0 mL
05109469

IVD 002 2-8 °C

LOT 00000000 2000-01-31

Pos.4 - 07598939700 - V2



Pos.5 - 07598955700 - V1



Pos.6 - 07598998700 - V1





Pos.7 - 07598971700 - V1



Pos.8 - 07599005700 - V1



IL6
CalSet
Elecsys
and **cobas e** analyzers
05109469 **590**
LOT 000000
 01



Pos.9 - 07599013700 - V1

Nhãn sản phẩm
IL-6 CalSet (05109469190)

Chất hiệu chuẩn xét nghiệm **định lượng** IL-6
(Interleukin-6)
IL-6 CalSet

Mã sản phẩm: 05109469190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:

Công ty TNHH Roche Việt Nam

Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,


Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , Số lô: xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991