

## *Document Detail*

---

**Type:** 2100-MS-S  
**Document No.:** 2100-MS-S200372[1]  
**Title:** Universal Tray Adesivo  
**Comment:**  
**Status:** CURRENT  
**Effective Date:** 19-Jun-2021

## *Approval*

<u>Owner Role</u>	<u>Sign-off By</u>	<u>Sign-off Date</u>	
2100-Document Owner Zhermack Export Document Owner	Stefania Quaglia	16-Jun-2021 1:58 pm	GMT
2100-R&D Director Zhermack R&D Director	Massimo Rossi	18-Jun-2021 12:08 pm	GMT
2100-QA/RA Manager Zhermack QA/RA Manager	Barbara Baruzzo	17-Jun-2021 8:50 am	GMT



### TRAY ADHESIVE SPECIFICATION

<b>ZHERMACK ITEM CODE:</b>	<b>C700025</b>	<b>PRODUCT NAME:</b>	<b>UNIVERSAL TRAY ADHESIVE</b>
<b>INTENDED USE:</b>	<b>ADHESIVE FOR METAL IMPRESS. TRAY</b>	<b>PRODUCT SHELF LIFE:</b>	<b>24 MONTHS</b>
<b>Test Conditions</b>	<b>Specifications</b>	<b>Standard Applied/Method</b>	
Room Temperature (°C)	( 23,0 ± 1 ) °C	2100-WI QC-001	
Relative Humidity (%)	( 50,0 ± 10 ) %	2100-WI QC-001	
Test Temperature (°C)	( 23,0 ± 1 ) °C	\	
<b>Features</b>	<b>Specifications<sup>1)</sup></b>	<b>Standard Applied/Method</b>	
Appearance	Viscous fluid-trasparent	2100-AM QC-004	
Colour	Cobalt Blue	2100-AM QC-004	
Traction resistance (3 tests average) N/mm <sup>2</sup>	> 50	2100-AM QC-038	

1) Values are referred to the features of the product up to the expiry date. The test is made in conformity to our internal methods: mixing ratio in weight, product hand mixed at humidity and temperature conditions mentioned in this specification. Standards and internal methods are the current version at the date of the test. Recommended STORAGE TEMPERATURE 5 ÷ 27°C and TERMS OF SHIPPMENT : -15 ÷ 40°C max 30 days and max 75%UR