

IFU – ULTRASOUND GEL

WARNING!

*Do not apply on injuries skin.
Do not use and get the
treatment if irritation occurs.*

Purpose for use

*External use only. This product is for
multiple use. As coupling medium
between transducer and external body
surface in ultrasound scanning
procedure.*

Instruction for use

*Apply a sufficient amount of gel on the
transducer or over the treatment area.
When treatment is completed,
remove the gel with towel or tissue
from the skin and device.*

Ingredients

*Purified Water, Carboxylic Polymer,
Triethanolamine, DMDM Hydantoin.*

Contraindication

*Not for internal use i.e. invasive
procedure or biopsy application.
Not for use with defibrillator.*