

# SENSA STAT CONTROL LEVEL - 3

## For use with

SCLBL1675

**Blood Gas Analyzers Only.**

## **Intended Use**

















Sensa Stat Controls (Level-2) used for monitoring of Sodium, Potassium, Ionized Calcium, pH, Chloride, Glucose, Lactate, pCO<sub>2</sub>, pO<sub>2</sub> on the Blood Gas Analyzer. Sensa Cores Sensa Stat Controls (Level-3) is not made for examine specimens from the human body for diagnostic purposes and made for clinical laboratory equipment use only. Sensa Cores Sensa StatControls (Level-3) is Single use only.

## Pack Size:

**Sensa Stat Level - 3 : 10 Ampules (10 x 1.8 mL)**

## **Storage and Stability**

Store at 18-25 °C until expiration date mentioned on labels

 CE	CE Marking	 IVD	In vitro diagnostic medical device
 Manufacturer	Manufacturer	 Caution	Caution
 Date of Manufacture	Date of Manufacture	 Do Not Use if Package is Damaged	Do Not Use if Package is Damaged
 Use by Date	Use by Date	 ISO	ISO Certified
 LOT	Batch Code	 Keep Dry	Keep Dry
 Keep Away from Sunlight	Keep Away from Sunlight	 Do not discard along with household waste	Do not discard along with household waste
 Storage & Operating Temperature	Storage & Operating Temperature	 Non Flammable	Non Flammable
 Consult instructions for use	Consult instructions for use	 Keep away from Flame	Keep away from Flame

***Irritant!*** Irritating to eyes and skin. Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

 **Sensa Core Medical Instrumentation Pvt. Ltd.**  
Plot No: 3, Export Promotion Industrial Park, IDA, Pashamylaram,  
Sangareddy (District) - 502307, Hyderabad, Telangana, INDIA.  
Tel: +91 8455 - 223400 www.sensacore.com

EC REP

Quality & Compliance  
SCLBL1675  
SCLBL1675  
SCLBL1675  
SCLBL1675

ISO 9001

ISO 13485