

Nhãn chính

MTTS CPAP CANNULA
SINGLE PATIENT USE ONLY

SIZE: 0
PART NO: CNL-10000

MTTS Co., Ltd
House No. 26, Alley 41, An Duong Vuong
Street, Tay Ho District, Hanoi City, Vietnam

2016-10
2021-10

LOT 16100009

REF CNL10000

DEHP-FREE
LATEX

Outside diameter (OD) selection:

size	OD/mm	code	approx weight of baby
0	3.1	•	1,500 grams and below
1	3.6	••	1,500 - 2,500 grams
2	4.4	•••	2,500 grams and above

Rev. 1.07 www.mtts-asia.com MADE IN MALAYSIA

Nhãn phụ

Gọng mũi đôi dùng cho máy trợ thở CPAP

Ký hiệu: MTTS Cannula

Sản phẩm dùng một lần.

Cỡ:

Tên và địa chỉ chủ sở hữu số Đăng ký lưu hành:

Công ty TNHH Chuyển giao Công nghệ và Dịch vụ Y tế
Số 26 ngõ 41 đường An Dương Vương, Phường Phú Thượng,
Quận Tây Hồ, TP Hà Nội.

Website: www.mtts-asia.com.vn +84 24 3766 6521

Sản xuất tại Malaysia

Ngày sản xuất: Xem trên nhãn chính

Hạn sử dụng: Xem trên nhãn chính

Số lô: xem LOT trên nhãn chính.

Thông tin về cỡ: Xem trên nhãn chính.

