

**MDK-5638K(IV) Electric Hospital Bed With Eight Functions** 电动医用病床(八功能)



**Technical parameters:**

Size: 2200\*1050\*500-950mm

**Standard Accessory:**

- Detachable ABS head/foot board with roller bumpers
- Full width : Including the side rails: 1050mm
- Bed base height range: Maximum height: 950mm & Minimum height: 500mm
- Patient touching surface: 1950\*900mm
- Back section panel: 0-75°
- Leg section panel: 0-45°
- Trendelenburg: 0-12°
- Anti-trendelenburg: 0-12°
- Left and right tilt: 0-12°
- Maximum safe load: 250KG
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- Embedded rails control
- ACP(Nurse controller)
- Center control lock
- Chassis with ABS cover

**Standard accessories:**

Nurse controller & Rail controller & Drainage & Detachable PP head and foot board & PP side rail

**Optional accessories:**

Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

**技术参数:**

•尺寸: 2200\*1050\*500-950mm

- 床高范围: 最高: 950mm 最低: 500mm
- 背板调节: 0-75°
- 腿部调节: 0-45°
- 整床向前倾斜: 0-12°
- 整床向后倾斜: 0-12°
- 整床向左倾斜: 0-12°
- 整床向右倾斜: 0-12°
- 最大承重: 250KG
- 可拆卸式ABS床头床尾板
- 可折叠式ABS护栏
- 护栏和床尾配备嵌入式控制器按钮
- 床体床面钢制磷化处理喷塑配备ABS软连接
- 可拆解式点滴架1支, 达型型尿袋构
- 最大承重: 250KG
- 四轮同步刹车系统1组
- 防尘式ABS底座

**标配:**

护士控制器、护栏控制器、引流挂构、可拆卸PP床头床尾板、PP护栏

**选配:**

床垫, 床边桌, 床头柜, 备用电池, CPR控制器



**MDK-5618K(II) Electric Hospital Bed With Five Functions** 电动医用病床(五功能)



**Specification:**

L2150\*W1050\*H500-800mm

**Standard Accessory:**

- Detachable ABS head/foot board
- Collapsible ABS guardrails
- Handle with lock
- Nurse control panel
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- Center control lock
- Chassis with ABS cover

**Optional accessory:**

Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

**规格:**

L2150\*W1050\*H500-800mm

**标配:**

- 可拆卸式ABS床头床尾板
- 可折叠式ABS护栏
- 可锁式控制器
- 挂置式护士控制器
- 床体床面钢制磷化处理喷塑
- 可拆解式点滴架1支, 达型型尿袋构
- 四轮同步刹车系统1组
- 防尘式ABS底座

**选配:**

床垫, 床边桌, 床头柜, 备用电池, CPR控制器



挂置式护士控制器  
Nurse control Panel



可锁式控制器  
Handle with lock