

MDK-5618K(III) Electric Hospital Bed With Five Functions 电动医用病床(五功能)



Specification:

L2150*W1050*H500-800mm

Standard Accessory:

- ▶ Detachable ABS head/foot board with roller bumpers
- ▶ Collapsible ABS guardrails
- ▶ Integrated stamping frame and mattress base are made of steel epoxy
- ▶ I.V. Pole & Drainage hook with PU cover
- ▶ Center control lock
- ▶ Chassis with ABS cover
- ▶ Handset 1pc

Optional accessory:

Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:

L2150*W1050*H500-800mm

标配:

- ▶ 可拆卸式ABS床头床尾板(带防撞轮)
- ▶ 可折叠式ABS护栏
- ▶ 床体床面钢制磷化处理喷塑
- ▶ 可调解式点滴架1支, 沾塑型尿袋钩
- ▶ 四轮同步刹车系统1组
- ▶ 防尘式ABS底座
- ▶ 手动控制器1个

选配:

床垫, 床边桌, 床头柜, 备用电池, CPR控制器

MDK-5618K(II) Electric Hospital Bed With Five Functions 电动医用病床(五功能)



Specification:

L2150*W1050*H460-750mm

Standard Accessory:

- ▶ Electric motor with manual crank together
- ▶ Detachable ABS head/foot board
- ▶ Collapsible ABS guardrails
- ▶ Integrated stamping frame and mattress base are made of steel epoxy
- ▶ Backrest with X-RAY examination tray
- ▶ I.V. Pole & Drainage hook with PU cover
- ▶ Center control lock
- ▶ Handset 1pc

Optional accessory:

Mattress & over-bed table & bedside cabinet & Back-up battery

规格:

L2150*W1050*H460-750mm

标配:

- ▶ 手电动组合控制系统(断电情况下正常使用)
- ▶ 可拆卸式ABS床头床尾板
- ▶ 可折叠式ABS护栏
- ▶ 床体床面钢制磷化处理喷塑
- ▶ 背部为X射线板功能
- ▶ 可调解式点滴架1支, 沾塑型尿袋钩
- ▶ 四轮同步刹车系统1组
- ▶ 手持控制器

选配:

床垫, 床边桌, 床头柜, 备用电池

MDK-5618K(IV) Electric Hospital Bed With Five Functions 电动医用病床(五功能)



Specification:

L2150*W1050*H460-720mm

Standard Accessory:

- ▶ Detachable ABS head/foot board
- ▶ Collapsible ABS guardrails
- ▶ Integrated stamping frame and mattress base are made of steel epoxy
- ▶ I.V. Pole & Drainage hook with PU cover
- ▶ Center control lock

Optional accessory:

Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:

L2150*W1050*H460-720mm

标配:

- ▶ 可拆卸式ABS床头床尾板
- ▶ 可折叠式ABS护栏
- ▶ 床体床面钢制磷化处理喷塑
- ▶ 可调解式点滴架1支, 沾塑型尿袋钩
- ▶ 四轮同步刹车系统1组
- ▶ 手持控制器

选配:

床垫, 床边桌, 床头柜, 备用电池, CPR控制