

MDK-5618K(V) Electric Hospital Bed With Five Functions 电动医用病床(五功能)



Specification:

L2150*W1050*H500-800mm

Standard Accessory:

- Detachable ABS head/foot board
- Embedded rails control
- Embedded nurse control on the foot board
- Collapsible ABS guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- Center control lock
- Chassis with ABS cover

Optional accessory:

Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:

L2150*W1050*H500-800mm

标配:

- 可拆卸式ABS床头床尾板
- 可折叠式ABS护栏
- 床尾板配备嵌入式控制器按键
- 护栏配备嵌入式控制器按键
- 床体床面钢制磷化处理喷塑
- 可调解式点滴架1支,达塑型尿袋构
- 四轮同步刹车系统1组
- 防尘式ABS底座

选配:

床垫, 床边桌, 床头柜, 备用电池, CPR控制器

Economic style

MDK-5611K(I) Electric Hospital Bed With Five Functions 电动医用病床(五功能)



Specification:

L2150*W950*H400-720mm

Standard Accessory:

- Detachable ABS head/foot board
- Collapsible aluminum alloy guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- Center control lock
- Hand remote control

Optional accessory:

Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:

L2150*W950*H400-720mm

标配:

- 可拆卸式ABS床头床尾板
- 可折叠式铝合金护栏
- 床体床面钢制磷化处理喷塑
- 可调解式点滴架1支,达塑型尿袋构
- 四轮同步刹车系统1组
- 手持控制器

选配:

床垫, 床边桌, 床头柜, 备用电池, CPR控制器

MDK-5611K(II) Electric Hospital Bed With Five Functions 电动医用病床(五功能)



Specification:

L2150*W950*H250-720mm

Standard Accessory:

- Detachable ABS head/foot board
- Collapsible steel coated guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- Center control lock
- Hand remote control

Optional accessory:

Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:

L2150*W950*H250-720mm

标配:

- 可拆卸式ABS床头床尾板
- 可折叠式铁制护栏
- 床体床面钢制磷化处理喷塑
- 可调解式点滴架1支,达塑型尿袋构
- 同步刹车系统1组
- 手持控制器

选配:

床垫, 床边桌, 床头柜, 备用电池, CPR控制器