

MDK-2611K Electric Hospital Bed With Two Functions 电动医用病床(二功能)



Specification:
L2150*W950*H500mm

Standard Accessory:

- Detachable ABS head/foot board
- Collapsible aluminum alloy guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- Center control lock
- Hand remote control

Optional accessory:

- Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:
L2150*W950*H500mm

标配:

- 可拆卸式ABS床头床尾板
- 可折叠式铝合金护栏
- 床体床面钢制磷化处理喷塑
- 可拆解式点滴架1支,沾塑型尿袋
- 四轮同步刹车系统1组
- 手持控制器

选配:

- 床垫,床边桌,床头柜,备用电池, CPR控制器

MDK-2611L(I) Electric Hospital Bed With Two Functions 电动医用病床(二功能)



Specification:
L2150*W950*H500mm

Standard Accessory:

- Detachable ABS head/foot board
- Collapsible aluminum alloy guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- 5寸装饰盖刹车轮
- Hand remote control

Optional accessory:

- Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:
L2150*W950*H500mm

标配:

- 可拆卸式ABS床头床尾板
- 可折叠式铝合金护栏
- 床体床面钢制磷化处理喷塑
- 可拆解式点滴架1支,沾塑型尿袋
- 5寸装饰盖刹车轮
- 手持控制器

选配:

- 床垫,床边桌,床头柜,备用电池, CPR控制器

MDK-2611L(II) Electric Hospital Bed With Two Functions 电动医用病床(二功能)



Specification:
L1900*W900*H500mm

Standard Accessory:

- Detachable ABS head/foot board
- Detachable steel coated guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V Pole & Drainage hook with PU cover
- Heavy duty caster with brakes 5"
- Hand remote control

Optional accessory:

- Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control

规格:
L1900*W900*H500mm

标配:

- 可拆卸式ABS床头床尾板
- 可拆卸式金属喷塑护栏
- 床体床面钢制磷化处理喷塑
- 可拆解式点滴架1支,沾塑型尿袋
- 5寸装饰盖刹车轮
- 手持控制器

选配:

- 床垫,床边桌,床头柜,备用电池

