

**MDK-2618L(I)** Electric Hospital Bed With Two Functions 电动医用病床(二功能)



**Specification:**

L2150\*W1050\*H500mm

**Standard Accessory:**

- Detachable ABS head/foot board
- Collapsible ABS guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V. Pole & Drainage hook with PU cover
- Heavy duty caster with brakes 5"
- Hand remote control

**Optional accessory:**

Mattress & over-bed table & bedside cabinet & Back-up battery

**规格:**

L2150\*W1050\*H500mm

**标配:**

- 可拆卸式ABS床头床尾板
- 可折叠式ABS护栏
- 床体床面钢制磷化处理喷塑
- 可调节式点滴架1支, 达型尿管架
- 5寸装饰盖刹车车轮
- 手持控制器

**选配:**

床垫, 床边桌, 床头柜, 备用电池

**MDK-1611L** Electric Hospital Bed With One Function 电动医用病床(一功能)



**Specification:**

L2150\*W950\*H500mm

**Standard Accessory:**

- Detachable ABS head/foot board
- Collapsible aluminum alloy guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V. Pole & Drainage hook with PU cover
- Heavy duty caster with brakes 5"
- Hand remote control

**Optional accessory:**

Mattress & over-bed table & bedside cabinet & Back-up battery

**规格:**

L2150\*W950\*H500mm

**标配:**

- 可拆卸式ABS床头床尾板
- 可折叠式铝合金护栏
- 床体床面钢制磷化处理喷塑
- 可调节式点滴架1支, 达型尿管架
- 5寸装饰盖刹车车轮
- 手持控制器

**选配:**

床垫, 床边桌, 床头柜, 备用电池

**MDK-2618L(II)** Electric Hospital Bed With Two Functions 电动医用病床(二功能)



**Specification:**

L2150\*W950\*H500mm

**Standard Accessory:**

- Detachable Head/Foot Board made of ABS plastic
- Collapsible ABS guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V. Pole & Drainage hook
- Heavy duty caster with brakes 5"

**Optional accessory:**

Mattress & over-bed table & bedside cabinet  
Spare battery & CPR control

**规格:**

L2150\*W950\*H500mm

**标配:**

- 可拆卸式ABS床头床尾板
- 可折叠式ABS护栏
- 床体床面钢制磷化处理喷塑
- 可调节式点滴架1支, 引流袋挂钩
- 5寸装饰盖刹车车轮

**选配:**

床垫, 床边桌, 床头柜, 备用电池, CPR控制器

**MDK-3011K(I)** Electric Hospital Bed With Three Functions 电动医用病床(三功能:家具型)



**Specification:**

L2100\*W980\*H400-720mm

**Standard Accessory:**

- Detachable solid wood head/foot board
- Collapsible aluminum alloy guardrails
- Detachable MDF guardrails
- Integrated stamping frame and mattress base are made of steel epoxy
- I.V. Pole & Drainage hook
- Center control lock
- Chassis with ABS cover
- Hand remote control

**Optional accessory:**

Mattress & over-bed table & bedside cabinet & Back-up battery

**规格:**

L2100\*W980\*H400-720mm

**标配:**

- 可拆卸式高密度床头床尾板
- 可折叠式铝合金护栏
- 可拆卸式高密度床边框
- 床体钢制磷化处理喷塑
- 可调节式点滴架1支, 达型尿管架
- 四轮同步刹车系统1组
- 防尘式ABS底座
- 手持控制器

**选配:**

床垫, 床边桌, 床头柜, 备用电池