

# Thermo Fisher Scientific (Asheville) LLC

Country of Origin: United States

Model Number: 88300V

Item Number: 132VM4C01L

Rev: 1.02

Suffix: 65

STG	REFRIGERANT			REFRIGERANT			OIL	
	TYPE	Kg	BAR	TYPE	Kg	BAR	TYPE	Kg
1	R-404a	0.595					POE	0.624
2	R-290	0.025	0.345	R-508b	0.312	9.515	AB	1.049
0								
0								

Volts: 230

Amps: 9.5

Hz: 50

Phase: 1

Test Pressure

Stage 1

High: 31.03 Bar

Low: 15.86 Bar

Stage 2

High: 27.58 Bar

Low: 15.86 Bar

Wt: 257 kg

Cap: 421 L

Climate: N

Temp: -66C

Serial Number



1160485501170124



This product is covered by one or more of the following US patents: 6408948, 6418131, 5702632, 6100672, 6119466, 6177271, 6216479, 6282914, 6297047, 6324893, 6324895, 6374620, 6374620, 6477489, 6481224, 6482637, 6543820, 6667356, 6631825. D494080 and/or pending and foreign patents.

Người đại diện hợp pháp của cơ sở



GIÁM ĐỐC

Vũ Tuấn Anh