



## Document Information:

Document Author: Schnabel, Helena {DOMC~MANNHEIM}  
Business Area / Unit: Roche Professional Diagnostics  
Confidentiality: Confidential  
Document Class: Device Master Record  
Document Type: Device Master Record  
Document Creator: Schnabel, Helena {DOMC~MANNHEIM}  
Document Lifecycle Status: Signed  
Valid From: 16-Jul-2024 13:55:54 (UTC)  
Valid To:  
Document Title: C\_0201\_RSP\_08058733000  
Document Number: 0000000000001200000302543  
Document Version: 05  
Template: No

Global Group: Regulatory Affairs  
Global SubGroup: Specification Prod. (Registr.)  
Local Group: Reg. Specification  
Site: RDG Germany  
Department: C\_DOMCC  
Document Description: Tina-quant Transferrin ver.2 (TRSF2) - "Only for use in  
Regulatory Affairs"  
Justification for change: Revision

## Electronic Signatures:

Signed By: schnabeh (Helena Schnabel {DOMCSD})  
Role: Reviewer  
Signature Differentiation: Quality Control  
Signed Date: 16-Jul-2024 13:03:22 (UTC)

Signed By: massothk (Katja Lienhard {DOMCC})  
Role: Approver  
Signature Differentiation: Quality Control  
Signed Date: 16-Jul-2024 13:21:45 (UTC)

# Specification

Cat. No.: 08058733190  
08058733214

Product: Tina-quant Transferrin ver.2 (TRSF2)

Parameter	Specification
-----------	---------------

## Reagents-working solutions

### Buffer solution (B):

pH value (25 °C):	7.00 - 7.40
-------------------	-------------

### Anti-human transferrin antibodies (C):

pH value (25 °C):	6.50 - 7.50
-------------------	-------------

## Performance test Serum application

### Recovery on cobas c 303 analyzer:

PreciControl ClinChem Multi 1	(Cat. No. e.g. 05947626190):	90 - 110 %
-------------------------------	------------------------------	------------

PreciControl ClinChem Multi 2	(Cat. No. e.g. 05947774190):	90 - 110 %
-------------------------------	------------------------------	------------

### Recovery on cobas c 503 analyzer:

PreciControl ClinChem Multi 1	(Cat. No. e.g. 05947626190):	90 - 110 %
-------------------------------	------------------------------	------------

PreciControl ClinChem Multi 2	(Cat. No. e.g. 05947774190):	90 - 110 %
-------------------------------	------------------------------	------------

**Parameter****Specification****Performance test Serum application****Recovery on cobas c 703 analyzer:**

PreciControl ClinChem Multi 1 (Cat. No. e.g. 05947626190): 90 - 110 %

PreciControl ClinChem Multi 2 (Cat. No. e.g. 05947774190): 90 - 110 %

**Performance test Urine application****Recovery on cobas c 303 analyzer:**

PreciControl ClinChem Multi 1 diluted (Cat. No. e.g. 05947626190): 85 - 115 %

PreciControl ClinChem Multi 2 diluted (Cat. No. e.g. 05947774190): 85 - 115 %

**Recovery on cobas c 503 analyzer:**

PreciControl ClinChem Multi 1 diluted (Cat. No. e.g. 05947626190): 85 - 115 %

PreciControl ClinChem Multi 2 diluted (Cat. No. e.g. 05947774190): 85 - 115 %

**Recovery on cobas c 703 analyzer:**

PreciControl ClinChem Multi 1 diluted (Cat. No. e.g. 05947626190): 85 - 115 %

PreciControl ClinChem Multi 2 diluted (Cat. No. e.g. 05947774190): 85 - 115 %

**Also required:**

Calibrator f.a.s. Proteins (Cat. No. 11355279216)

**Roche Diagnostics GmbH  
Quality Control Mannheim**Sandhofer Strasse 116  
D-68305 Mannheim  
Phone +49-621-759-0