



Product description – Vacuukit

The Vacuukit is intended to mix and deliver bone cements for cementoplasty & arthroplasty. It is composed of a reservoir, a funnel, a mixer plunger, a tube and an injection handset.

The Vacuukit has been designed to contribute to the porosity reduction during the preparation of bone cement by applying vacuum and thereby increase the mechanical stability.

Intended use

The Vacuukit is intended to mix and deliver bone cements for cementoplasty & arthroplasty.

Composition

Components	Material
Stand	ABS
Funnel	Polypropylen
Mixer	Polycarbonate, polypropylene, polyethylene, silicon
Injection tube	LDPE
Injection handset	Polycarbonate + GF, Stainless steel, ABS
Vacuum tube with filter	PVC, polypropylene, Charcoal

Reference

Designation	Reference number
Vacuukit	T060416



Product Description – S5

The S5 system is intended to mix and deliver PMMA bone cements. It is composed of a shaker with a metal mixing ball, a funnel, a cement delivery gun and an extension tube.

S5 kit is packaged in a box and sealed in a double barrier packaging which guarantees its sterility. The blister is sterile and can be put on the operating table. Single use only.

Indications

The S5 system is intended to mix and deliver bone cements for vertebroplasty, kyphoplasty and cementoplasty.

Composition

- ✓ Injector gun: Polycarbonate, AISI 316, AISI 302, Aluminium, ABS, EPDM, POM, Polyamide.
- ✓ Extension tube : Polyamide, AISI 316, EPDM, AISI 316L
- ✓ Shaker : PP, EPDM
- ✓ Funnel : ABS
- ✓ Ball : Glass

References

Designation	Reference number
S5	T060410