

SHANGHAI SHENWANG RADIATION PROTECTIVE EQUIPMENT CO., LTD.

Date: Mar 13, 2017

Technical Standard

DESCRIPTION	TESTING ITME	STANDARD VALUE	RESULT
X-radiation Protective Glass Plates ZF2 All Sizes	Lead Equivalent	1.76 mmPb	OK
	Transparent Rate	85%	OK
	Specific Gravity	4.12 g/cm ³	OK
	Stability	A	OK
	Peak Voltage	120 kV (Min. 1.64 mmPb) 150 kV (Min. 1.53 mmPb)	OK
	Physical Characteristics	Conforming to Standard	OK

Legal Representative: Gu Di

Signature:



SHANGHAI SHENWANG RADIATION PROTECTIVE EQUIPMENT CO., LTD.