

## CA32

### Fundamental Parameters

Dimension	Diameter: 13.1cm; Height: 17.1cm
Weight	135.6g

### Composition and Material

Cover lid	PP
Canister body	PS
Yellow Adapter	PE

### Package

25PCS/CTN	477mm*462mm*347mm	6.3kg
-----------	-------------------	-------

### Performance Parameters

Volume	800ML
Vacuum resistance	80kPa
Service Life (est.)	2 years

### Storage Conditions

Ambient Temperature	-10°C~50°C
Relative Humidity	5%~85% R.H
Storage Pressure	700hPa~1060hPa

### Operating Conditions

Operating Temperature	0°C~40°C
Relative Humidity	30%~75% R.H
Atmospheric Pressure	700hPa~1060hPa

### Accessory Specification (Optional)\*

Solidifier	Absorbent Polymer Resin
------------	-------------------------

**\*Note:**

All the accessories are for single use only.

Please use the recommended product or contact us.

