

# MẪU NHÃN GIƯỜNG ELEGANZA SMART JUNIOR



Address of Manufacturer  
**CÔNG TY CỔ PHẦN**  
**TRANG THIẾT BỊ VÀ**  
**CÔNG TRÌNH Y TẾ**  
 M.M.D.N. 0102778235 - C.T.C.P.  
 LINE 1 spol. s r.o.  
 Zelebrice 5, 274 01 Štávy  
 Czech Republic  
<http://www.linet.com>

1D Bar code EAN128  
 (Serial Number)

Manufacturing Date  
 (Year-Month-Day)

DI (Device Identifier)  
 GTIN (Global Trade  
 Item Number)

2D Bar Code  
 (GS1 Data  
 Matrix)  
 DI+PI=UDI

MODEL: 1XXYYYYY-kYY (Name of Bed)  
 SERIAL NO.: 20160063476  
 ORDER NO.: xx/yyyy/zzzz YYYY-MM-DD  
 (01)00592654120521  
 (11)150906  
 (21)20160063476  
 =250 kg =165 kg =210 kg V<sub>n</sub>: 230 V, 50/60 Hz  
 P: 370 VA (0,37 kW) IP X4  
 INT.: 10%, Max. 2min/15 min

Symbols

Serial Number

Electrical Specification

PI (Product Identifier)