

Bassinet Cart (PB-1100, PB-1200)



Photograph is PB-1100

- Clear plastic bassinet
Dimensions (mm): (755 L x 438 W x 250 H)
- Nameplate holder
- 75 mm dia. caster with brakes on opposite corners
- Polyurethane foam mattress
Dimensions (mm): (650 L x 350 W x 50 TH)
- Wire basket with powder coating finishing



Colour	
Pink	Blue
PB-1100	PB-1200



Pediatric Bed (PB-2100, PB-2200)



Photograph is PB-2100

- Crank-operated back raise adjustment
- Steel sliding side rails
- Pressed steel bottom
- 100 mm dia. Paramount original caster with brakes on opposite corners
- Polyurethane foam mattress
Dimensions (mm) : (1660 L x 830 W x 100 TH)
- I.V.Hanger Rod (Option) can be attached



Back Raise
0° - 75°

Colour	
Pink	Blue
PB-2100	PB-2200



Psychiatric Bed (PB-3101, PB-3102, PB-3103, PB-3104)



Photograph is PB-3102

- 7 pressed steel bottoms are welded to 3 long bars
- Each bottoms have holes for restraint belt
- Bed frame is supported by four legs and additional pipe leg in the middle
- Head and foot board are fixed with screws
(Board type can be chosen)

Board Type			
Steel Pipe		Plastic	
Without Corner Bumper	With Corner Bumper	Green	Beige
PB-3101	PB-3102	PB-3103	PB-3104

