



# SERIM MEDICAL EQUIPMENT CO., LTD.

118, Daerim-ro, Nami-myeon, Seowon-gu, Cheongju-si, Chungcheongbuk-do, Republic of Korea

E-mail: serimkorea1@gmail.com Mobile. 82-(0)10-6413-5487

## SPECIFICATION

MODEL	SR-DC8550-2S	NAME	DRESSING CART	ORIGIN	SOUTH KOREA
-------	--------------	------	---------------	--------	-------------



**Round Bumper**

※The design of the product(Caster) may change without notice in order to improve quality.

### Technical Specifications

1. Overall Size(W x D x H): 850 x 500 x 900mm
2. Material: Made of stainless steel304
3. Two shelves with **4side guard rail**
4. 4" Swivel casters (Swivel-2ea, Brake 2ea)

### Features

1. Serim Dressing Cart has an open design that is easy to clean and dry, load and maneuver for containing all supplies necessary for changing dressings of surgical or injured patients.
2. It is constructed of stainless steel 304 for rust proof and non contamination.
3. Tabletop edges are safety hemmed.

### Consist of

- 1.Main body : 1Set