

AG-XD205 Blood donor chair



Seat width: 600mm

Seat width, including armrests: 880mm

Total length: 1880mm

Total weight: 70kg

Height: 530mm

Caster: dia.100mm

Back cushion: 190*490*120mm

Cushion: PU foam

Upholstery: PVC leather

Frame: Q235 steel

Safe maximum load: 240kg

Back segment angle: -12° -75°

Trendelenburg: -12°

Leg segment angle: -80° -12°

AEGEAN'S AFTER SALE



As long as purchasing AEGEAN products,from the date on finish and certify products

Warranty:18 months

Warranty for motors:36 months

After-sale service:15 years

In the warranty period:

Due to non-human factors,the damaged products will receive our free spare parts with technical guidance.

Beyond the warranty period

We charge the cost of spare parts but provide technical guidance and other after-sale service for free

Whole-hearted service everyday!

Final interpretation belongs to AEGEAN

Tel: 0086-512-56382825 Fax:0086-512-58170725 E-mail:info@aegeanchina.com aegean@aegeantech.com

website : www.aegeantech.com www.aegeanchina.com www.aegeachina.com

Aoyang International Building,Zhangjiagang,Jiangsu,215600,China

JIANGSU AEGEAN TECHNOLOGY CO.,LTD