

# SMACO

AN EMS COMPANY



Swiss Technology



**NEW**

## SYG-200 Amalgamator

- MCU Control - accurate and reliable.
- Special capsule clamping device.
- Working time can be adjusted freely (01-99 seconds).
- Can be used for blending amalgam and glass ionomer capsules.

### Specifications

- Power input: AC220V/50Hz
- Power consumption:  $\leq 220V A$
- Mixing frequency: 4300T/min
- Timer: 01-99 seconds
- Noise level:  $\leq 65 dB$
- Weight: 4.0kg
- Dimension: 250 x 220 x 145mm

Product Code: SMAC-SYG200

**1 Year  
Warranty!**



Since the founding of SMACO in the 1990's, SMACO has made it's name in the market as one of the leading manufacturers of curing lights and amalgamators.