



Date: June 13, 2018

Dear Valued customer,

Product : Latex Examination Glove Powdered , White, Textured at finger tip, Non-Sterile**Certificate of Analysis**

Characteristics	Specification (TH-LXMD-1FB)	Average Actual result Latex Powdered (TH-LXMD-1FB)
1.Length (mm.)		
- XS	Min 220	241
- S	Min 220	241
- M	Min 230	241
- L	Min 230	242
- XL	Min 230	-
2. Width (mm.)		
- XS	< /=80	78
- S	80+/-10	83
- M	95+/-10	94
- L	110+/-10	107
- XL	> /=110	-
3.Single Wall Thickness (mm)		
- Finger	Min 0.11	0.12
- Palm	Min 0.08	0.10
4. Physical Properties		
4.1 Before Aging		
- Elongation (%)	Min 650	776
- Force (N)	Min 7.0	9.0
4.2 After Aging		
- Elongation (%)	Min 500	743
- Force (N)	Min 6.0	8.6
5. Powder residue on powder gloves	6-15 mg/dm ²	9.3
6. Protein	Max 150 ug/g	118
7. Hole	AQL < /= 2.5	

Remark: Average test results in January-December 2017

Brand Name	Specification
SATORY(20X100)(STA)-LXM-01FB	TH-LXMD-1FB

Best Regards,

Ms. Vanlinee Laohachaiyakul
Product Manager