



**OL-400**

**OPERATION LAMP(X7CM)**



**SPECIFICATIONS**

- 7등
- Housing Diameter 780mm
- Illumination Intensity 1~6 Level Adjustment 150,000Lux
- Lamp LED
- Bulb 12V-168W
- Power Consumption 170W
- Total Weight 45kg
- Minimum Ceiling Height 2,500mm

**OL-300**

**OPERATION LAMP(X5CM)**



**SPECIFICATIONS**

- 5등
- Housing Diameter 560mm
- Illumination Intensity 1단~6단 조절가능 140,000Lux
- Lamp LED
- Bulb 12V-120W
- Power Consumption 120W
- Power Source AC 220V
- Total Weight 41kg

**OL-500 OPERATION LAMP(5+5CM,DUOS)**



**SPECIFICATIONS**

• Total weight 66kg

**5등(Sub)**

- Maximum Light Output 140,000Lux
- Housing Diameter 560mm
- Color Temperature(K) 4,300
- Color Rendering Index(CRI) > 92
- Light Field Diameter(cm) 22 ~ 28
- Number of LED Bulbs 40
- Average Bulb Life(hours) 40,000
- Minimum Requirement for Ceiling Height 2,300mm
- Power Consumption 120W
- Power Input AC 220V 50/60Hz

**5등(Sub)**

- Maximum Light Output 140,000Lux
- Housing Diameter 560mm
- Color Temperature(K) 4,300
- Color Rendering Index(CRI) > 92
- Light Field Diameter(cm) 22 ~ 28
- Number of LED Bulbs 40
- Average Bulb Life(hours) 40,000
- Minimum Requirement for Ceiling Height 2,300mm
- Power Consumption 120W
- Power Input AC 220V 50/60Hz

**LED BULB LIFE SPAN**

• The advantage of new LED light with character of high performance and lower consumption, its bulb service life over 40000 hrs is the most efficient and economical compared with halogen bulb.

**INTRODUCTION**

Focus & light field adjustable (Aluminum sterile hand)

- Adjustment of the focus and diameter of light field by sterilizable handle(Aluminum)

**POWER LED CONTROL**



- On/Off Switch
- Electronic light intensity control