



PMMA Disk Vita16/w0.5(bleaching)/Clear

For temporary crowns and bridges



- Open Type
Diameter: 98.5mm (80~120mm available)
Thickness: 8~25mm
- Zirkonzahn Type
Diameter: 95mm (80~120mm available)
Thickness: 8~25mm
- Amann Girrbach Type
Diameter: 98mm (80~120mm available)
Thickness: 8~25mm

Wax Disk

Wax Disk (Color Green) for wax patterns to casting

Wax Disk Alpha (Color Grey) for wax patterns to rapid casting

- Open Type
Diameter: 98.5mm (80~120mm available)
Thickness: 8~25mm
- Zirkonzahn Type
Diameter: 95mm (80~120mm available)
Thickness: 8~25mm
- Amann Girrbach Type
Diameter: 98mm (80~120mm available)
Thickness: 8~25mm




PMMA Block A1,A2,A3

For temporary crowns and bridges

- PMMA Block for Roland DG DWX-4 (76×40×20mm)
- PMMA Block Pin Type
- PMMA Block 2 Pin Type

Wax Block

- Wax Block for Roland DG DWX-4 (76×40×20mm)




CAD/CAM Milling Bar Line-up 2 types

For DWX-50 and DWX-4 of Roland DG

Diamond Coated Type

-10 times hardened durability compared to titanium coated bars



Diamond Coated Type/DLC Type


Type	Shape	Under head Length(mm)	Shank Diameter(mm)	Overall Length(mm)
RADIUS Ø 3mm	RADIUS	16	4	45
Ø 2mm	BALL	16	4	45
Ø 0.8mm	BALL	12	4	45
Ø 0.4mm	BALL	6	4	45

DLC (Diamond Like Coated) Type

-2 times hardened durability compared to titanium coated bars



Diamond Like Coated Type- **approx 2 times**, Diamond Coated Type- **approx 10 times** more durable compared to general coated bars



Increase a processing efficiency by using a RADIUS end milling bar

By utilizing a RADIUS end milling bar, a larger area can be milled that results in the increase of processing efficiency and combined utilization of both RADIUS end milling bar and Ball end milling bar Ø2mm enables you to shorten the processing time by 30% compared to the time of single application of Ball end milling bar Ø2mm.



YAMAHACHI DENTAL MFG. CO.
54-1 OCHIGARA NISHIURA-CHO, GAMAGORI-CITY, AICHI PREF. JAPAN

Tel: (+81) 533-57-7121 / Fax: (+81) 533-57-1764

E-mail: box@yamahachi-dental.co.jp

Web: www.yamahachi-dental.co.jp

