

New Products EXAMINATION TABLE series

FLAT TYPE


1985 (L) x 750 (W) x 650 (H) mm



Photograph is PK-1102

PK - 1102	Caster Type
PK - 1105	Leg Type

- Flat Table Type
- Two Table Tops with Polyurethane Foam
Mattress covered with vinyl leather
- Steel Frame (Powder Coated)
- Fixed Steel Pipe Board (Chrome Plated)
- Option : (1) **TENTE** Precision Casters
(2) Leg Type with Rubber Cap



BACK RAISE TYPE

1985 (L) x 750 (W) x 650 (H) mm




Photograph is PK-1112



PK - 1112	Caster Type
PK - 1115	Leg Type

- Back raise adjustment
- Two Table Tops with Polyurethane Foam
Mattress covered with vinyl leather
- Steel Frame (Powder Coated)
- Fixed Steel Pipe Board (Chrome Plated)
- Option : (1) **TENTE** Precision Casters
(2) Leg Type with Rubber Cap

Back raise angle can be controlled by means of a yellow gas spring lever (0 ~ 90°)




BACK RAISE AND HI-LO TYPE

1985 (L) x 750 (W) x 550 ~ 870 (H) mm




Photograph is PK-1122



PK - 1122	Caster Type
PK - 1125	Leg Type

- Back raise and hi-lo adjustment
- Two Table Tops with Polyurethane Foam
Mattress covered with vinyl leather
- Steel Frame (Powder Coated)
- Fixed Steel Pipe Board (Chrome Plated)
- Option : (1) **TENTE** Precision Casters
(2) Leg Type with Rubber Cap

Hi-Lo handle



OPTION

CASTER TYPE

With TENTE Precision Casters (Corrosil)



Double ball bearing swivel caster enhances mobility and easy maneuverability. Zinchromatized and additionally palted Corrosil housing enables exceptional durability against rust corrosion.

Ø 125 mm Casters (2 Casters with brakes on diagonal)

LEG TYPE

Leg with Rubber Cap




For not moval usage

I.V. POLE

KC - 56

I.V. Hanger Rod


- I.V. Hanger Rod (Stainless Pipe Rod) with attachment part
- 1030 ~ 1775 (H) mm



PC - 5003

I. V. Stand

- I.V. Stand with Casters (40mm dia)
- 437(W) x 460(D) x 1290 ~ 2075(H) mm



SAFETY BELT

Safety belts are available subject to special order

PILLOW

Pillow is available subject to special order

PT. PARAMOUNT BED INDONESIA



MM2100 INDUSTRIAL TOWN, BLOCK M-1-1
EXPORT PROCESSING ZONE, CIKARANG BARAT
BEKASI 17520, INDONESIA
Tel.: 62-21-8980715, Fax.: 62-21-8980716
URL <http://www.paramount.co.id>
Email : marketing@paramount.co.id



PARAMOUNT BED CO.,LTD.

OVERSEAS TRADE DIVISION
14-5, 2 CHOME HIGASHISUNA, KOTO-KU, TOKYO
136-8670 JAPAN

TEL: 81-3-3648-2961 FAX: 81-3-3648-2420
URL <http://www.paramount.co.jp>
Email : overseas@paramount.co.jp

Effective from November 3, 2006